

## SPECIAL PACKAGE OFFER

### *Urology procedures are safer than ever*

The Great White<sup>®</sup> Platinum Stirrups provide safe & quick positioning assistance for improved surgical site access while offering exceptionally precise placement.

#### **Features**

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -30° low to 75° high; Adduction -9° and abduction +25°
- 600 lb. (272 kg) surgical patient weight capacity
- The raised lateral fin boots effectively minimize the possibility of nerve damage
- Included Platinum Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs



### Special Platinum Stirrup Package #8194

- SchureMed Platinum Stirrups
- Uro-Trapper System
- SpringLoc Rail Clamp
- Stirrup Transport & Storage Dolly

#### **All Products Also Available A La Carte**

- SchureMed Platinum Stirrups
  - SpringLoc Clamp
  - Uro-Trapper System
  - Uro-Trapper Drain Bags, 10/cs
  - Stirrup Transport & Storage Dolly
  - Schure Socket Clamp, US
  - Universal Rail Clamp, US
  - Disp. Clear View Leg Drapes, Sterile, 10/cs\*
- \*Sterile products are non-returnable.*

