

# 100 SERIES

## *Procedure Chair*

- Powered height, backrest & tilting
- Open base design for great ergonomics
- 550 lb patient capacity



Shown in Taupe fabric with Black Aesthetic Upgrade, Arm Rests and Adjustable Head Rest Platform & Pad.



FM584901

made in the USA  
with US & imported parts



OAKWORKS<sup>®</sup>  
**MEDICAL**  
*Position for Success*

# 100 SERIES *Procedure Chair*

The 100 Series Procedure Chair gives you the flexibility to perform many different procedures while supporting great ergonomics and patient comfort. The twin tower design supports up to 550 lb and the top is designed for supine or prone positioning.

Shown in standard black powder coated base and Seafoam fabric



## 100 SERIES PROCEDURE CHAIR #67395-T-E\* 2.5" Comfort #73406-T-E\* 4" Comfort

FEATURES	
3 Section Top	27" (69cm) wide x 74" (188cm) long with Head Rest [80" (203cm) long with Adjustable Head Rest extended]
Back rest	27" (69cm) wide x 26" (66cm) long
Seat	27" (69cm) wide x 24" (61cm) long
Leg rest	22" (56cm) wide x 13" (33cm) long
Motion 1	Electric Height Range: 19" - 35" (48-89cm)
Motion 2	Electric Trendelenburg/Reverse Tilt: +15°/-5°
Motion 3	Electric Backrest: 0-80°
Motion 4	Manual 4 position Leg Rest 0-90°
Foot Control	Adjust Height, Trendelenburg/Reverse & Backrest
Weight Capacity	550 lb (250kg)
Product Weight	300 lb (136kg)
Padding	2.5" (6cm) or 4" (10cm) Comfort Foam™
Paper Roll Holder	Accommodates 18" (46cm) roll (roll not included)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate

OPTIONS <i>must be ordered with table.</i>		#
Aesthetic Upgrade	Acrylic base cover and skirt (specify color)	73177
Arm Rests	Fold down design	73174
Hand Control	Controls all functions (can be used with foot control)	69928
Non-USA electric system	Please contact us for order details	73134
Procedure Tray	12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters) Available only with 4" Comfort Foam™ Reduces foam thickness in seat section to 2"	73042
Stirrups (Retractable Foot Rests)	Pull out from under seat section, fully adjustable	73175
T-rail, Fixed	Available for factory mount on one or both sides of the seat. 14" (36cm) long	73043
ACCESSORIES <i>order any time. See accessories on page 59.</i>		#
Adjustable Head Rest Platform	Adjustable angle & height, removable	71899
Procedure Chair Face Rest Crescent	For use when prone positioning is required. Use with Adjustable Head Rest Platform	69876-T*
Head Rest Pad	10" (25cm) long x 10" (25cm) wide Use with Adjustable Head Rest Platform	71976-T*
Large Head Rest	Adjustable in length, removable 11.5" (29cm) long x 20" (51cm) wide	72717-T*
T-rail Adapter	Attaches to chair at multiple locations for use with USA format T-rail attachable accessories.	61902

\*Fabric code (T) and Enclosure color code (E) (if ordering Aesthetic upgrade) must be specified when ordering.

## ADJUSTABLE HEAD REST PLATFORM

#71899

Make your patient comfortable with this head rest platform that provides adjustable angle & height to fit their needs.



## HEAD REST PAD

#71976-T\*

Use the head rest pad with the adjustable head rest platform for optimum comfort.

10" (25cm) long x 10" (25cm) wide



ACCESSORY

## LARGE HEAD REST #72717-T\*

This accessory is available by itself without the adjustable head rest platform and the head rest pad. It is removable and is adjustable in length.


11.5" (29cm) long x 20" (51cm) wide



ACCESSORY

## 100 SERIES ARM REST #73174

Fold down arm rests help the patient feel more secure. They fold down completely out of the way to ensure great access.



OPTIONAL

## HAND CONTROL #69928

Add a hand control to adjust Trendelenburg, height, and backrest.



OPTIONAL

## STIRRUPS #73175

Adjust to a wide range of angles & lengths to accommodate any patient or procedure. Retractable design stores inside the table.



OPTIONAL

## PROCEDURE TRAY #73042


Retractable, removable stainless steel. 12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters)



OPTIONAL

## T-RAIL ADAPTER #61902

Multiple locations allow you to choose the best spot for this easily installed & removed accessory capable of handling up to 150 lb (68kg). Can be used with any American standard accessory. 8" (20cm) long




ACCESSORY

### ENCLOSURE COLORS


Choice of 4 acrylic colors to match your decor. \*Color code must be specified.



### PADDING




2.5" Comfort Foam™  
Comfort and support



4" Comfort Foam™  
Superior comfort and support

### TERRATOUGH™ FABRIC

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.



We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.

CONTACT:



*View our full  
line of Medical  
products*



#72703  
Scan this code  
with your  
smart phone

717.235.6807 [www.oakworksmed.com](http://www.oakworksmed.com) email: [customerservice@oakworks.com](mailto:customerservice@oakworks.com)