

PF SERIES

Exam Table

- Maximum comfort for extended treatments
- Beautiful wood base with choice of finishes
- Up to 4 powered motions



PF400
shown in Seafoam fabric
with optional Ebony
Lacquer and Locking
Casters



PF200



PF250




PF400



ISO
13485
Medical Devices
Quality
Management

FM584901

made in the USA
with US & imported parts



OAKWORKS[®]
MEDICAL
Position for Success

PF SERIES

Exam Table

INCREDIBLE COMFORT FOR LONGER TREATMENTS



Electric 4 Section Top



MODEL PF400

- #70882-T* 31" W 4" Comfort
- #73309-T* 29" W 4" Comfort

- Motion 1 Electric Height: 22"-38" (56-97cm)
- Motion 2 Electric Tilt: 20°/20°
- Motion 3 Electric Backrest: 0-53°
- Motion 4 Electric Leg rest



Lift-Assist Backrest Top



MODEL PF250

- #70889-T* 31" W 3" Comfort
- #73409-T* 31" W 4" Comfort
- #73454-T* 29" W 3" Comfort
- #73566-T* 29" W 4" Comfort
- #73583-T* 27" W 3" Comfort
- #73584-T* 27" W 4" Comfort

- Motion 1 Electric Height: 22"-38" (56-97cm)
- Motion 2 Electric Tilt: 20°/20°
- Motion 3 Lift Assist Backrest: 0-80°

Flat Top



MODEL PF200

- #70890-T* 31" W 3" Comfort
- #73410-T* 31" W 4" Comfort
- #73567-T* 29" W 3" Comfort
- #73568-T* 29" W 4" Comfort
- #73585-T* 27" W 3" Comfort
- #73586-T* 27" W 4" Comfort

- Motion 1 Electric Height: 21"-37" (53-94cm)
- Motion 2 Electric Tilt: 20°/20°

PF SERIES EXAM TABLES

FEATURES	
Top Dimensions	See above for widths x 73" (185cm) long
Foot Control	Controls all powered movements
Weight Capacity	550 lb (250kg)
Product Weight	235 lb (106kg)
Padding - PF400	4" (10cm) Comfort Foam™
Padding - PF250 & PF200	3" (8cm) or 4" (10cm) Comfort Foam™
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS <i>must be ordered with table.</i>		#
4" (10cm) Locking Casters	To easily move your table in the room	69926
**Wood Finishes - see pg. 78	Choice of 4 stains or 4 lacquer finishes	73136
Non-USA electric system	Please contact us for order details	73134
ACCESSORIES <i>order any time.</i>		#
Arm Rests (PF400)	Auto-adjusting & removable for easy access	74734-T*
Universal Table Extender with 4" Padding	21" (53cm) W x 12.5" (32cm) L	74736-T*
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support.	PKG7554*
Paper Roll Holder	Vanilla - Accommodates 21" (53cm) roll (roll not included) Not FDA listed	11580

*Wood finish and fabric color code must be specified.

The PF Series Exam Tables are designed to maximize comfort for extended treatments while offering the beauty of wood to create a Spa-like atmosphere. The fully electric design comes with a super soft and supportive padding system.

CASTERS
#69926



OPTIONAL

HAND CONTROL

Add a hand control to control all functions.



enquire for part number

OPTIONAL

ADJUSTABLE SIDE ARM RESTS #66760-T*

These incredibly versatile, removable arm rests rotate horizontally and always stay level. 7" (18cm) wide at widest point x 23.5" (60cm) long



(PF 400 Only)

ACCESSORY

PAPER ROLL HOLDER

#11580 vanilla #11349 black

This simple design is adjustable to accommodate different size rolls up to 21" (53cm) wide & plugs into the table's outlet holes. (paper roll not included)



ACCESSORY

QUICKLOCK™ HEAD REST & AERO·CEL PAD

#PKG4186-T*

The platform allows precise cervical flexion and height positioning over a wide range for both prone and supine treatments.



ACCESSORY

UNIVERSAL TABLE EXTENDER


#68697-T* 3" padding
#68735-T* 4" padding

Perfect as a head rest or foot rest, it has two rods that fit into outlets on the table. Extends table to 84". 21" (53cm) W x 12.5" (32cm) L




ACCESSORY

PADDING



3" Comfort Foam™
More comfort and support
(Not available on PF 400)



4" Comfort Foam™
Superior comfort and support

STAINS & LACQUERS

Oakworks[®] finishes are highly water resistant and can be wiped down with disinfectants as needed. Stains allow the natural grain of the wood to show through and create depth. Lacquers are more opaque and have a more solid appearance.



Optional Stains




Optional Lacquers



TERRATOUGH™ FABRIC

Meets the requirements of ISO 10993-1 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.



We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.

CONTACT:

717.235.6807 www.oakworksmed.com email: customerservice@oakworks.com