

Benefits

FOR THE SONOGRAPHER

- Build-A-Bed concept enables you to customize the table with all the features you need
- Design supports good body mechanics, minimizing strain and injury.
- Unobstructive wheel base design allows mobility around the table & prevents injury

FOR THE PATIENT

- Ultimate in comfort with our 3” padding
- Easy on & off, wheelchair accessible
- Top design provides security for all shapes & sizes, with a 550 lb. bariatric rated lift & load capacity

FOR THE PURCHASER

- Customizing your table is easy & cost effective, only pay for the features you want
- Low maintenance, all electric design
- FDA listed to ensure safety & compliance
- Industry leading warranty
3 years parts, 2 years labor

*For product details,
please refer to our
Medical Catalog*



Scan this code
with your
smart phone



OAKWORKS[®]
MEDICAL

*Things to consider
when choosing a table:*

- What types of ultrasound procedures will be performed currently and in the future?
- What electric movements are a must for a new table?
- Will the table be moved around the room or mainly positioned in one spot?
- What accessories will be required?
- Are there any ergonomic issues to be addressed with the purchase of a new table?
- What is the budget?

*Contact us for help finding the
best solution for you!*

Contact:

717.235.6807 www.oakworksmed.com

71778 rev 08.25.15


OAKWORKS[®]
MEDICAL
Position for Success

Helping you select the right table!

ULTRASOUND

- Industry Leading Warranty
- Ergonomically advanced
- Customizable
- Bariatric strength

made in the USA



Excellent ergonomics reduce sonographer and patient injury.

Perfect for...

- Heart & Vascular Labs
- Hospitals/Medical Centers
- Women's Health Centers
- Private Practice

Build-A-Bed

"Build-A-Bed" concept lets you customize Oakworks exam tables to the individual needs of each ultrasound suite. You only pay for the features you want & need. It's an innovative way of ordering a truly custom built ultrasound exam table.

Start with any of our General Ultrasound tables, then simply add the features you want based on the type of ultrasound specialty you perform. Or you can order one of our Specialty tables that is preconfigured to meet your needs.

Build-A-Bed™ is Trademarked by Sound Ergonomics.

OAKWORKS[®] MEDICAL Position for Success

GENERAL SPECIFICATIONS

- Height Ranges from 22" - 38" (56-97cm)
- Weight Capacity: 550 lb. (250kg)
- Choose from 27" or 30" (69-76cm) Widths
- Hand Control Standard, Optional Foot Control
- Locking Casters
- Paper Roll Holder
- **Industry Leading Warranty**

OPTIONS & ACCESSORIES (vary by model)

- Patient Lift Compatible Base
- Central Base Locking System
- Retractable Foot Rests (*stirrups*)
- Adjustable Access Back Support
- Sliding Cardio Panel
- Side Rails
- Procedure Tray
- Accessory Mounting Plate
- Detachable Arm Board
- T-Rail Adapter
- Table Extender
- Face Rest

TERRATOUGH™ FABRIC CHOICES




*Inspired by the
leading ergonomic
experts in sonography*

made in the USA



ISO 13485:2003
FM 584901



VASCULAR
WITHOUT FOWLER



VASCULAR
WITH FOWLER



GENERAL
RECTANGULAR TOP



GENERAL
WITH FOWLER



GENERAL
3-SECTION



WOMEN'S IMAGING



MULTI-SPECIALTY



ECHOCARDIOGRAPHY